

THERAPEUTIC CARDIAC RHYTHM MONITORING



POWERHEART[®] CRM[™]
cardiac rhythm module

ACCESSORIES



ESTABLISHING A NEW STANDARD OF CARE FOR CARDIAC PATIENTS WITH THE POWERHEART CRM AND ACCESSORIES

Introducing the Cardiac Science Powerheart CRM Line of Accessories.

Featuring the Powerheart CRM Transport Carrying Case specifically designed to facilitate the use of the CRM during transport when you need to "Pick Up & Go."

- Water-Repellant Exterior
- Detachable adjustable shoulder strap with cushioned shoulder pad
- 2 side pockets accommodate & provide easy access to electrodes, power supply, cables & Quick Reference Guide



The unique design of the Powerheart CRM Transport Carrying Case allows the CRM to be continuously plugged into AC without having to remove the CRM from the carrying case.

Other accessories for the Powerheart CRM include the Wall (GCX-compatible) and IV Pole Mounting Brackets which provides additional CRM placement options to fit a variety of space configurations and requirements.

To learn more about Cardiac Science's unique Powerheart CRM and its line of accessories, visit our website www.cardiacscience.com, or email info@cardiacscience.com, or call 1.866.289.5649 (USA and Canada only).



CARDIAC SCIENCE

CREATING HEARTSAFE ENVIRONMENTS

ACCESSORIES

ELECTRODES

Powerheart CRM Disposable Therapy Electrodes

P/N: 170-2007

- Reliable ECG monitoring, defibrillation, and external pacing performance
- Specifically designed for use with the Powerheart CRM Cardiac Rhythm Module
- Integrated reference electrode offers enhanced discrimination from motion artifact, reducing potential for false alarms during patient movement
- Compact design and radiotranslucent material allows for clearer views of the chest under X-ray and fluoroscopy
- Packaged 10 sets per box



ECG Monitoring Electrodes (3-pack)

P/N: 170-0124

- Standard ECG monitoring snap-on electrodes
- For use with Powerheart CRM ECG Cable
- Packaged 80 pouches per case



CABLE

Powerheart CRM Patient Therapy Cable P/N: 162-0058-003

- 8 ft. (2.44 m)
- For use with Powerheart CRM Disposable Therapy Electrodes
- Compatible only with the Powerheart CRM

Powerheart CRM ECG Cable (3-Lead) P/N: 162-0057

- AHA labeled version
- For use with standard ECG monitoring electrodes
- Allows additional 3-lead ECG monitoring through the Powerheart CRM

Powerheart CRM ECG Pass-Thru Cables

P/N: 170-2088-001

- Allows output of ECG signal to compatible Philips®/Agilent patient monitoring systems

P/N: 170-2089-001

- Allows output of ECG signal to compatible GE® Marquette patient monitoring systems

TRANSPORT SYSTEM

Powerheart CRM Transport Carry Case

P/N: 168-1020-001



For more information contact Cardiac Science at:

Corporate Headquarters

1900 Main Street, Suite 700, Irvine, CA, 92614, USA
Phone: +1.949.797.3800 Fax: +1.949.951.7315

Minnetonka Operations

5474 Feltl Road, Minnetonka, MN, 55343, USA
Phone: +1.952.939.4181 Fax: +1.952.939.4191

International Operations

Kirke Vaerloesevej 14, DK-3500 Vaerloese, Denmark
Phone: +45.4438.0500 Fax: +45.4438.0501

Cardiac Science Holdings (UK) Limited

Unit 13, Crown Industrial Estate, Canal Road, Timperley
Altrincham, Cheshire, WA14 1TF, United Kingdom
Phone: +44.161.926.0000 Fax: +44.161.926.0002

Cardiac Science Deutschland GmbH

Max Planck Strasse 7B
D-52249 Eschweiler, Germany
Phone: +49.2403.5046.53 Fax: +49.2403.5046.54

Worldwide Order Entry:

Toll Free: +1.800.991.5465 (USA and Canada)
Phone: +1.949.797.3806
Fax: +1.949.851.4418
Email: aedorders@cardiacscience.com

USA/Canada Customer Service/Technical Support:

Toll Free: +1.888.466.8686
Phone: +1.949.797.3807
Fax: +1.949.851.4425
Email: aedsupport@cardiacscience.com

International Customer Service/Technical Support:

Phone: +44.161.926.0000
Fax: +44.161.926.0002
Email: aed@cardiacscience.co.uk

MOUNTING SYSTEMS

Powerheart CRM Wall-Mount System P/N: 168-2003

- Rotating support arm extension mount provides secure mounting of Powerheart CRM
- Compatible with most hospital wall-mount systems

Powerheart CRM IV Pole Mount System P/N: 168-2011

- Clamp design provides secure mounting of Powerheart CRM
- Compatible with most hospital I.V. poles

LITERATURE

Powerheart CRM User/Service Manual Binder, English P/N: 112-2000-001

- Complete operating and maintenance instructions for Powerheart CRM

Powerheart CRM Quick Reference Guide, English P/N: 112-2028-001

- Short instructions for use for the Powerheart CRM

MISCELLANEOUS

Powerheart CRM ECG Training Simulator

P/N: 168-2001

- ECG/Pacing simulator designed for use with Powerheart CRM
- For demonstration and training purposes only



Powerheart CRM Defibrillation Test Verifier

P/N: 168-2005

- Incorporates display light that flashes during Powerheart CRM energy discharge

Powerheart CRM Event Review Kit, English P/N: 109-2000-001

- Includes Powerheart CRM-to-PC service port cable and software (CD-ROM)
- Allows for downloading, storing, and reviewing a patient's cardiac events

Powerheart CRM AC Power Adapter P/N: 139-0186

- AC power supply for use with the Powerheart CRM

Powerheart CRM Power Cord P/N: 167-0102-001

- Shielded power cord compatible for North America

Powerheart CRM Replacement Battery P/N: 170-2077

- Sealed lead acid battery pack

Powerheart CRM Extended Warranty Service Contract

- Covers product for an additional 18 months (months 19-36)



CARDIAC SCIENCE

Nasdaq: DFIB

1.866.289.5649

WWW.CARDIACSCIENCE.COM

EMAIL: INFO@CARDIACSCIENCE.COM



Specifications subject to change without notice. CRM, FirstSave, IntelliSense, MasterTrak, MDLink, Powerheart, RescueLink, RescueReady, RHYTHMx, STAR, and Survivalink are trademarks and registered trademarks of Cardiac Science, Inc. Microsoft and Windows are registered trademarks of Microsoft Corporation. All other trademarks are the property of their respective owners.